


*Equipment:*

# **ELECTROCARDIOGRAPH**



**Cardio Técnica**  
Argentine Manufacturer of Hospital Equipment




**Model: RG-401plus**

## Equipment:

# ELECTROCARDIOGRAPH

Model: RG-401plus



## Accessories

- Carrying bag
- Patient cable
- Power cord
- Ground wire connection
- Kit with 4 electrode holder clips
- Kit with 6 precordial electrodes
- Gel bottle
- Thermal paper roll
- User's Handbook

## Technical Specifications

Model: RG-401plus

- Mode Selection: 1-3 channels
- Digital membrane keyboard
- High resolution thermal/dot-matrix printer
- Two modes of operation: automatic and/or manual
- LCD screen for data displaying and of the cardiac signal
- PC connection interface. Software optional
- Two printing speeds: 25 or 50 mm/sec
- Printing mode selection
- Sensitivity: 0.5, 1 or 2 cm/mV. Auto gain
- Optically insulated power source
- Cabinet made from molded injected plastics of high mechanical rigidity
- Power requirements: 110/220 V
- Rechargeable batteries (charger included)
- Three-language menu
- Recording paper: 50mm x 30m C
- Calibration: 1mV +/-5%
- Filters: 50/60Hz
- Leads (12): I-II-III-AVR-AVL-AVF-V1-V2-V3-V4-V5-V6
- Patient circuit: Preamplifier optocoupled with floating input; defibrillation discharge protection
- Dimensions: 10.63" x 7.87" x 3.15" (270 x 200 x 80 mm)
- Weight: 3.91 lbs. (1,80 kg) (batteries included)

