


Equipment:


# PULSE OXIMETERS



**Cardio Técnica**  
Argentine Manufacturer of Hospital Equipment



**Model: OXI-401**



**Model: OXI-411**



Equipment:

# PULSE OXIMETERS

Model: OXI-401

## Technical Specifications

Model: OXI-401

- LED digital display.
- Rechargeable batteries included.
- Programmable, visual and audible alarms.
- Small, lightweight and portable body.
- Wide variety of sensors: adult, pediatric and neonatal.


- Digital membrane keyboard.
- Parameters: oxygen and pulse saturation.
- Data obtained from oximeter sensors.
- Floating patient input.
- No connection to ground or power.
- Oxygen saturation indicator: 0% to 100%
- Heart rate indicator: 0% to 100%
- Alarm setting:

Upper limit: 0 to 250 ppm

Lower limit: 0 to 250 ppm

Upper limit: 100%

Lower limit: 50%







Equipment:

# PULSE OXIMETERS

Model: OXI-411

## Technical Specifications

Model: OXI-411



- On-screen plethysmographic waveform.
- Digital LED display.
- Rechargeable batteries (12 hours of continuous operation).
- Programmable, visual and audible alarms.
- Small, lightweight and portable body.
- Wide variety of sensors: adult, pediatric and neonatal.
- Alarm readings: 3-digit display per parameter.
- Alarm limit variation: linear, digital.
- QRS display: adjustable acoustic beep volume.
- Power supply: 110/220 V~
- Batteries: 12 V / 800 mAh
- Dimensions: 8.66" x 7.87" x 3.54" (220 x 200 x 90 mm).
- Weight: 2.42 lbs (1.10 kg).
- Mode of representation: On-screen continuous plethysmographic waveform with freezing capability & highly luminescent LED display (OXI-411 & OXI-401).

Product approved by Ministerio de Salud y Ambiente A.N.M.A.T. PM - 1230 - 2